




## Automated Fluorescence Immuno Assay System

# AFIAS-3

- Automated POC system
- Highly sensitive result based on advanced fluorescence technology
- Easy-to-use cartridge system
- LIS/HIS compatible



- AFIAS-3 is a fluorescence immunoassay-based POC analyzer that is an integrated Sample-to-Answer system featuring 3 parallel tests.
- AFIAS-3 is a compact POC system that provides fast and reliable result on the spot.
- AFIAS All-in-one cartridge enables convenience to users.



Specialized for  
panel tests in lab



ER or ICU where  
quick result is required



Convenience with  
All-in-One cartridge

## All-in-One cartridge for AFIAS platform

- All the reagents for assay contained
- Over 50 test parameters
- Ready-to-use cartridge



## Specifications

Model	AFIAS-3
Dimensions (L x W x H)	365 x 276 x 314 mm
Weight	9.3 kg
Analyzer input	100~240V AC, 50~60Hz 1.8A
User Interface	7" Touch-screen LCD
Data Output	Internal Printer/RS-232 port
Operating temperature	15°C to 35°C
Connectivity	RS-232 3 USB port / LAN port SD card slot

Automated Fluorescence Immuno Assay System

# AFIAS-3



\* AFIAS-3 will be released in Q1 of 2022.