






ichroma™ COVID-19 Ag test

Specific COVID-19 Antigen test


- Fast and precise result (12 mins)
- Small POC device
- Room temperature storage
- Nasopharyngeal swab or nasal aspirate specimens
- Connect to LIS/HIS
- Objective results






ichroma™ COVID-19 Ag test


• Test procedure



• High sensitivity (Time resolved fluorescence) immunoassay



• Clinical performance



* Based on the clinical data with ichroma™ II Influenza A+B reagent

• Specification

Product name	ichroma™ COVID-19 Ag
Assay type	TRFLFA*, COVID-19 Antigen
Sample type	Nasopharyngeal swab, nasal aspirate specimens
Sample prep / extraction	Collecting with cotton swab
Target	COVID-19 Antigen
Platform	ichroma™ II Reader
Unit	25 kit/box

• Ordering Information

Product	Cat. No.
ichroma™ COVID-19 Ag	CFPC-115
ichroma™ II Reader	FPRR021

* Note. A cotton swab is excluded in the components of the cartridge box.