



# ichromα™ COVID-19/Flu Ag Combo

## Testing for COVID-19 and Influenza A/B

- Small POC or hand-held device
- Room temperature storage
- Nasopharyngeal swab
- Connect to LIS/HIS
- Objective results






# ichroma™ COVID-19/Flu Ag Combo

- Health officials caution about the **twindemic** with COVID-19 and Influenza in coming flu season.
- **Both COVID-19 and influenza** can have common symptoms with various ranges from asymmetric to severe symptoms.
- **ichroma™ COVID-19/Flu Ag Combo** can be diagnosed for distinguishing of COVID-19 and Influenza in patients with symptoms.


## • Test procedure



## • Specification

Product name	ichroma™ COVID-19/Flu Ag Combo
Assay type	TRFLFA*
Sample type	Nasopharyngeal swab
Sample prep / extraction	Collecting with nylon swab
Target	COVID-19 Ag / Influenza A Ag / Influenza B Ag
Platform	ichroma™ II Reader / ichroma™ M2 Reader
Unit	25 kit / box

\* TRFLFA: Time-resolved fluorescent lateral flow assay



ichroma™ COVID-19/Flu Ag Combo

## • Ordering Information

Product	Cat. No.
ichroma™ COVID-19 Ag / Flu Ag Combo	CFPC-117
ichroma™ II Reader	FPRR021
ichroma™ M2 Reader	FPRR031

\* Note. A nylon swab is excluded from the components of the cartridge box.

