

When your heart needs...

AFIAS/ichroma™ BNP

- Fast processing
- Accurate result
- Simple procedure



For diagnosis and assessment of severity of heart failure

B-type natriuretic peptide (BNP) is secreted primarily by the ventricular myocardium in response to wall stress, including volume expansion and pressure overload.


• Specification



Item	ichroma BNP	AFIAS BNP
Sample type	WB / Plasma	
Anti-coagulant	EDTA	
Sample volume	35 µL	100 µL (30 µL, C-tip)
Time	12 min	
Working range	5 - 5,000 pg/mL	
Reference value	100 pg/mL	
Storage	20 month, 4-30 °C	
Catalogue No.	CFPC-121	SMFP-83

• Performance

iCHROMA-II



AFIAS-6

